

Lumena™ (Bendable-Tip) Ureteral Access Sheath

Provides reliable access to the ureter while protecting surrounding tissues and enhancing procedural control.



360° Lateral Flexibility

All-around flexibility for seamless maneuverability


270° Bending Capability

Maximum bending angle of 270° for precise control and access

**Highly compatible with
flexible endoscopes**

**Enhanced maneuverability
and precise control**

Lumena™ (Bendable-Tip) Ureteral Access Sheath



Product Specification

Model Specification	Effective Length L (mm)	Outer Diameter OD (mm)	Inner Diameter ID (mm)
UAS951135U	350	3.85	3.20
UAS951140U	400		
UAS951145U	450		
UAS111335U	350	4.45	3.85
UAS111340U	400		
UAS111345U	450		
UAS121435U	350	4.70	4.05
UAS121440U	400		
UAS121445U	450		
UAS131535U	350	5.00	4.40
UAS131540U	400		
UAS131545U	450		

*Full model specifications, options and reordering details can be found on our website or by reaching out to our support team.

Customizable Items

☐ Inner/outer packaging and accessories

☐ Product appearance and batch number

☐ Product specifications
(outer diameter, length, etc.)

☐ Catheter structure

☐ Bending length and angle

☐ Others



CatheGenix Pte. Ltd.

+65 64655909

info@CatheGenix.com

www.cathegenix.com

2 Tukang Innovation Grove,
#08-09/10 JTC Medtech Hub, Singapore 618305

